



## NJ1 & FREAK compact bikes

### THE NJ1 & FREAK KOMPAKTBIKES – RIDING FUN FOR BIG & SMALL

The NJ1 compact bike and the FREAK compact bike are all about riding enjoyment, flexibility and customisation. Their modular design allows the drive unit and chassis to be quickly detached – ideal for transport and storage. Precise steering geometry, a long wheelbase and balanced weight distribution ensure excellent handling characteristics. The NJ1 is aimed at ambitious leisure cyclists, whilst the FREAK is specifically tailored to the needs of children and teenagers. Both models can be individually configured to suit everything from touring to sport cycling.



#### #MaximumFlexibility

The NJ1 & FREAK compact bikes can be dismantled into two parts without tools: Drive unit and chassis. You can separate them in a few simple steps and reduce the pack size.



#### #CustomisedFitting

With an extensive range of options and accessories, your NJ1 or FREAK compact bike can be perfectly tailored to your needs. Suitable options are also available for users with limited mobility.



#### #ShiftAsYouLike

You can choose between hub or derailleur gears, as well as mechanical and electronic gear systems – operated by hand or with your chin.



#### #ChildrenAndYouthHandbike

FREAK compact bike: Light, compact and easy to handle – perfectly tailored to the needs of young riders.



NJ1



FREAK



## NJ1 & FREAK – AS INDIVIDUAL AS YOU ARE!

Your **NJ1** or **FREAK** compact bike is customised to your individual needs. You decide on a wide range of options – from the bottom bracket support, shifting and braking systems to the backrest, neck support and grip variants. Whether it's a touring bike or a sports bike: your NJ1/FREAK adapts perfectly to you. Your handbike, your choice.

### OPTIONS FROM PRO ACTIV:



Modular design: The drive unit and chassis can be separated in just a few seconds without tools, making transport and storage easier.



Individual adjustment options for the foot-rest, backrest and headrest ensure optimum comfort and ergonomics.



Reliable brakes, maximum control. With disc brakes or the optional coaster brake, you're in safe hands with both models.





Perfect gearing. Both models offer a wide selection of gears, from derailleur gears and hub gears to bottom-bracket gears. Whether mechanical or electric – you decide!

## TECHNICAL SPECIFICATIONS

|                                  |                                                                                            |
|----------------------------------|--------------------------------------------------------------------------------------------|
| Special features:                | Stable with excellent handling, with a range of specifications from touring to sports bike |
| Total weight:                    | NJ1 compact bike: from 15.8 kg<br>FREAK compact bike: from 14.8 kg                         |
| Pack size NJ1 compact bike:      | Chassis from 104 x 43 x 39 cm<br>Drive unit from 105 x 43 x 66 cm                          |
| Pack size FREAK compact bike:    | Chassis from 102 x 43 x 37 cm<br>Drive unit from 99 x 43 x 50.8 cm                         |
| Drive wheel & impeller sizes:    | 20"   26"                                                                                  |
| Frame lengths:                   | 102 cm   116 cm                                                                            |
| Seat widths:                     | 35 cm   39 cm   43 cm                                                                      |
| Max. standard load weight:       | up to 100 kg payload                                                                       |
| Max. load weight for Heavy Duty: | individual weight limit up to 120 kg<br>[NJ1 compact bike only]                            |

Data subject to design changes.

## FOLLOW US ON SOCIAL MEDIA



**PRO** **ACTIV**<sup>®</sup>

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11  
D-72359 Dotternhausen

info@proactiv-gmbh.de

[www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)

Rev. 08/2026\_01

Subject to technical modifications, printing errors & colour deviations.